


How to use an Automated External Defibrillator (AED)




- 1 Turn on the AED
- 2 Listen to the instructions
- 3 Attach the pads to the chest
- 4 Remove your hands from the patient before completing the next step
- 5 Press 'shock' only if the AED says 'shock advised'
- 6 Continue chest compressions until an ambulance arrives



Warning:

Make sure you and others stay clear when delivering the 'shock.'



نحن نهتم
We Care

مرکز جونز هوبكنز
أرامكو الطبي
Johns Hopkins
Aramco Healthcare